

NPS-XY-100A & NPS-XY-100A-UHV

100 μ m x 100 μ m Super Invar Stage

Nanopositioning Stage

The NPS-XY-100A and the NPS-XY-100A-UHV NanoMechanisms are piezo scanned flexure guided stage with integrated capacitance position sensors capable of sub-nanometer resolution and reproducibility. The NPS-XY-100A-UHV has been designed to operate down to 10⁻⁹ Torr.

Finite element analysis of the flexure guidance mechanisms has reduced parasitic angular motions to less than 25 micro-radians over the full 100 μ m range.

The Super Invar construction (CTE 0.3ppm K-1, compared to Aluminium at 23ppm K-1) minimizes thermal drift, which on a nanometer scale can be very important.

The unique iso-static mounting system ensures that stresses from the mounting system are properly relieved and establishes the center of the stage as the co-ordinate reference point.



Key features

- >100 μ m travel in each axis with sub-nanometer resolution
- Typically <0.005% hysteresis and <0.01% linearity error
- First resonant frequency >350Hz
- In-situ scanning and stepping response optimization
- Robust and reliable
- Super Invar construction

Typical applications

- • AFM, SPM, NSOM
- • High Precision Microscopy

Suggested controllers

- NanoScan NPC-D-6330 Series Controllers

Delivering Low noise, Low drift, High power, High Resolution and High position update rates.

Designed specifically to control Queensgate's Nanometer Precision Mechanisms incorporating capacitive sensors. They give precise positional feedback delivering highresolution and linearity of movement. The fast update rate and Queensgate control algorithms contribute to high speed positioning accuracy for dynamic applications that require high speed movement of the stage.

The PC software facilitates user optimisation of all operating parameters, including PID and notch filter set up. There are eight programmable slots, three which are populated to provide fast, medium and slow PID settings, the addition five slots are available for application specific settings.

Function playback provides user defined pre-programmed waveforms for applications such as raster scanning or constant velocity scanning. The calibration and dynamic settings are held in the stage eeprom which allows controllers to be interchanged with minimal performance changes.

Technical specification NPS-XY-100A

Parameter	Symbol	Value			Units	Comments
Static physical						
Material		Super Invar (Bright nickel plated)				
Size		100 x 100 x 23			mm	Note 1
		Minimum	Typical	Maximum		
*Range closed loop	d _{xp} -max	± 50	± 55		µm	
*Scale factor error (1σ)	δ _{bx1}			0.1	%	
Static stiffness			1		N·µm ⁻¹	
Resonant frequency:	0g load	F0-0	350		Hz	
	50g load	F0-50	260		Hz	
	1000g load	F0-100	120		Hz	
Maximum load				1	Kg	Note 2
Dynamic physical (Typical values)						
		Fast	Medium	Slow		Note 3
*3dB Bandwidth	B _{x-p}	53	20	4		Typical
*Small signal settle time	t _{xs-s}	15	30	130	ms	Note 4
*Position noise (1σ)	δ _{xp-n}	0.5	0.4	0.25	nmRms	Note 5
Error terms						
		Minimum	Typical	Maximum		
*Hysteresis (peak to peak)	δ _{xp-hyst}		0.005	0.015	%	Note 6
*Linearity error (peak)	δ _{xp-lin}		0.005	0.01	%	Note 6
*Rotational error	δφ _x		10	25	µradians	Note 7
*Rotational error	δθ _x		5	10	µradians	Note 7
*Rotational error	δγ _x		5	10	µradians	Note 7
Orthogonality	δθ _{orth}		8	10	mradians	

Notes

*These parameters are measured and supplied with each mechanism

- With 40mm diameter central aperture.
- Depends on orientation. 1 kg is the maximum load for gravity acting in the Z direction. 0.5 kg is the maximum load for gravity acting in the X or Y axes. Loads greater than 5 kg can cause damage to the flexure mechanism.
- For dynamic operation the servo-loop parameters are preset for different performances; the parameters are user settable via software control. Fast means the fastest the stage can stably move with less than 50 grams load. Medium means the maximum stable speed for loads up to 200 grams. Slow means the speed at which the servo loop is stable for all masses up to the maximum allowed mass – equivalent to low noise setting.
- Step and Settle time is the time taken to settle to within 2% of the step measured using an interferometer. The step settle time is a function of the servo loop parameters which are user controllable. The test step size is 2000 nm.
- The actual position noise of the stage. Measured with a laser interferometer sampling 1 Hz to 25 kHz.
- Percent of the displacement.
- Percent error over the full range of motion

General Note: Specifications assume standard 2m cable length, for custom cable lengths and UHV position noise, linearity and hysteresis may vary.

UNITED KINGDOM

Prior Scientific Instruments Ltd.
Units 3-4 Fielding Industrial Estate
Wilbraham Road, Fulbourn
Cambridge, CB21 5ET
United Kingdom
Email: inquiries@prior.com
Phone: +44 (0)1223 881711

U.S.A.

Prior Scientific, Inc.
80 Reservoir Park Drive
Rockland, MA. 02370
U.S.A.
Email: info@prior.com
Phone: +1 781.878.8442

GERMANY

Prior Scientific Instruments GmbH
Maria-Pawlowna-Str. 4
D-07743, Jena, Germany
Email: jena@prior.com
Phone: +49 (0) 3641 24 20 10

JAPAN

Kayabacho 3rd Nagaoka Bldg 10F,
2-7-10, Nihonbashi Kayabacho, Chuo-Ku,
Tokyo103-0025, Japan
Email: info-japan@prior.com
Phone: 03-5652-8831

CHINA

Prior Scientific Instruments (Suzhou) Ltd.
Room 1812, Honghai Building,
72 Xingdu Street, Suzhou Industrial Park,
Suzhou, 215000 China
Email: info-china@prior.com
Phone: +86 (0)512 6617 5866

