


USED ON

IF IN DOUBT ASK

DO NOT SCALE



PART MUST BE FREE FROM BURRS


DRG. NO. OEM-0030-I

2M STANDARD CABLE TERMINATES IN A  
15 WAY "D" TYPE CONNECTOR.  
MIN BEND RADIUS :-14mm

EMC CAP

TOP VIEW

DOWEL HOLE  
Ø2.007<sup>+0.005</sup><sub>-0</sub> × 5.0 DEEP  
IN 2 POSITIONS




UNDERSIDE VIEW

DOWEL HOLE  
Ø2.007<sup>+0.005</sup><sub>-0</sub> × 5.0 DEEP  
IN 2 POSITIONS

CABLE EXIT

HOLE Ø3.5 THRO'  
C'BORE Ø6.5 × 5.0 DEEP  
4 POS'NS  
FIXING TORQUE 60dNm TO 70dNm

HOLE M2.5 × 5.0 DEEP  
4 POS'NS  
FIXING TORQUE 40dNm TO 50dNm



5	ECR-0064-R	JR	JR						
4	NPS-7256-R	JR	CF	JR	PMR	MJM	13/04/99	ALL DIMENSIONS IN MM	
3	OEM-7004-R	YX/GJ	ASM	GMJ	PMR	PMR	29/10/98	TOLERANCES UNLESS OTHERWISE STATED	
2	OEM-7001-R	YX/GJ	GMJ	GMJ	PMR	MJM	10/02/98		
1		YX/GJ	ASM	GMJ	MJM	MJM	26/11/97	1 DECIMAL PLACE ±0.25 2 DECIMAL PLACES ±0.15 TRUE POSITIONAL ZONE DIA. 0.1 ANGULAR (DEGREES) ±0.5	
ISSUE	ECR No	DESN	DRAWN	CHKD	REA	APP'D	DATE		

© COPYRIGHT 1997-2004  
QUEENSGATE INST. LTD.

**Queensgate Instruments Limited.**  
WOODLAND ROAD, TORQUAY, DEVON, TQ2 7AY.



MATERIAL  
STAINLESS STEEL  
NON MAGNETIC

TITLE  
NPS-X-15A  
(INSTALLATION DRAWING)

FINISH  
CLEAN

SCALE	2-1	DRG. SIZE	A2	ISSUE	5
DRAWING NUMBER	OEM-0030-I			SHEET OF	1 / 1