


HOLE  
 M2.5 x 4.0 DEEP  
 EQUI-SPACED ON A 25.00 PCD  
 3 POSN'S  
 MAX. FIXING  
 TORQUE = 40/50cNm

HOLE  
 Ø1.60 x 3.00 DEEP  
 SUITABLE FOR DOWELS  
 2 POSN'S

HOLE  
 M1.6 x 3.0 DEEP  
 EQUI-SPACED ON 15.00 PCD  
 3 POSN'S



DO NOT CLAMP ON SLEEVE.  
 ONLY MOUNT THIS DEVICE  
 USING THE TOP OR BOTTOM  
 INTERFACE PLANES.

2M CABLE  
 2.80mm O/D  
 MIN. BEND RAD. :-10mm

15 WAY "D"TYPE PLUG

HOLE  
 M2.5 x 6.0 DEEP  
 EQUI-SPACED ON A 29.00 PCD  
 IN 3 POSN'S.  
 MAX. FIXING  
 TORQUE = 40/50cNm

HOLE  
 Ø1.60 x 3.00 DEEP  
 SUITABLE FOR DOWELS  
 2 POSN'S

Ø10.00 CLEAR APERTURE

|       |            |       |       |       |     |       |          |
|-------|------------|-------|-------|-------|-----|-------|----------|
| 2     | NPS-7193-R | YX/JR | JR    | JR    | PMR | MJM   | 23/2/99  |
| 1     |            | YX/JR | JR    | JR    | MJM | PDA   | 02/02/99 |
| ISSUE | ECR No     | DESN  | DRAWN | CHKD. | REA | APP'D | DATE     |
| 1     |            | 2     |       |       |     | 3     |          |

ALL DIMENSIONS IN MM  
 TOLERANCES  
 UNLESS OTHERWISE STATED

|                           |       |
|---------------------------|-------|
| 1 DECIMAL PLACE           | ±0.25 |
| 2 DECIMAL PLACES          | ±0.15 |
| TRUE POSITIONAL ZONE DIA. | 0.1   |
| ANGULAR (DEGREES)         | ±0.5  |

|   |   |              |                 |
|---|---|--------------|-----------------|
| © COPYRIGHT 1997<br>QUEENSGATE INST. LTD. | QUEENSGATE INSTRUMENTS LIMITED<br>WATERSIDE PARK, BRACKNELL, BERKSHIRE, RG12 1RB. |              |                 |
| MATERIAL                                  | TITLE INSTALLATION DRAWING:-<br>NPS-Z-15H   |              |                 |
| FINISH                                    | SCALE 1 : 1   | DRG. SIZE A3 | ISSUE 2         |
|   | DRAWING NUMBER NPS-0039-I   |              | SHEET 1<br>OF 1 |
| 1   | 2   | 3            | 8               |