







QGSP301XR - SP MULTIWELL HOLDER XR  
HOLDS ANSI SLAS STANDARD MICROPLATE (127.76 x 85.48 mm) (HOLDER MASS ~210g)



NOW INCLUDES PLASTIC 4 SLIDE INSERT  
HOLDER INSERT PART NO W4416





QGSP302XR - SP UNIVERSAL HOLDER XR  
HOLDS SINGLE 1"x3" (25x75mm) OR DOUBLE SLIDES 2" x 3" (50 x 75mm)  
OR DISHES 1"-3" (25 TO 75mm) DIAMETER (HOLDER MASS ~ 155g)






MOUNT DISH HOLDERS INSERT  
(SOLD SEPARATLY)  
FOR ENHANCED DISH MOUNTING  
AVAILABLE IN A RANGE OF SIZES




PART NO:  
QG-UNIV-DISH-35 (33-35mm DISH)  
QG-UNIV-DISH-40 (36-40mm DISH)  
QG-UNIV-DISH-55 (50-55mm DISH)  
QG-UNIV-DISH-60 (55-60mm DISH)




QGSP303XR - SP SLIDE HOLDER 1"x3" (25x75mm) & 2"x3" (50x75mm) HOLDER XR  
HOLDS EITHER SINGLE SIZE OR DOUBLE WIDTH SLIDES (HOLDER MASS ~165g)



QGSP304XR - FOUR SLIDE HOLDER XR  
HOLDS FOUR STANDARD 1"x3" (25mmx75mm) SLIDES (HOLDER MASS ~170g)



TITLE		SAMPLE HOLDER INFO SP STAGE SERIES				APRVD.		3RD ANGLE		DO NOT SCALE IF IN DOUBT ASK
DRG. No.	INSTALLATIONS NOTES		SHEET 1 OF 1	ISS.	A	DRAWN	J.C.			
MATL.	ALUMINIUM BODY				SCALE	1:1	TOLERANCE UNLESS OTHERWISE STATED ±0.1mm			
FINISH	BLACK ANODISE				DATE	22.05.25				
PRIOR Scientific Instruments Ltd. Fulbourn, Cambridge.										

ECN	REV	CREATED	APP	DATE	DETAILS
1st Rel	A	J.C.			1st Release