





CABLE 2m STANDARD

4 x Ø 3.30 ± 15
M4 - 6H ± 12
✓ Ø 4.20 X 90°, Near Side
MAX SCREW TORQUE 0.6N.m

4 x Ø 5.40 THRU



AXIS OF MOTION
(IF MOUNTED VERTICALLY
THIS SHOULD BE UPWARDS)

TITLE	INSTALLATION OF NPS-Z-250B / NPS-X-250B PIEZO STAGE	APRVD.	
DRG. No.	QG-0534-IN-100002	SHEET 1 OF 1	ISS. A
MATL.	ALUMINIUM STAGE BODY	DRAWN	J.C.
FINISH	NICKEL	SCALE	1:2
PRIOR Scientific Instruments Ltd. Fulbourn, Cambridge.		DATE	14.05.25
		3RD ANGLE	DO NOT SCALE IF IN DOUBT ASK
		TOLERANCE UNLESS OTHERWISE STATED	±0.1mm
		BREAK ALL SHARP EDGES	

ECN	REV	CREATED	APP	DATE	DETAILS
1st Rel	A	J.C.			1st Release