



NOTES  
WHEN USING THE STAGE LIFTING AGAINST GRAVITY THE STAGE MUST BE MOUNTED WITH THE ARROW  
UPWARDS TO MAINTAIN PRESSURE ON THE PIEZO ACTUATOR WITHIN.

TITLE	XP200 200um	APRVD.		3RD ANGLE		DO NOT SCALE IF IN DOUBT ASK
DRG. No.	QGXP200	DRAWN	J.C.	SCALE	2:1	TOLERANCE UNLESS OTHERWISE STATED ±0.1mm
MATL.		DATE	03.02.25			BREAK ALL SHARP EDGES
FINISH						
PRIOR Scientific Instruments Ltd. Fulbourn, Cambridge.						

ECN	REV	CREATED	APP	DATE	DETAILS
1st Rel	A	J.C.			1st Release